




# Kappy Pic Solution

**CE**  
**MD** 1936

Illustrations solely for the purpose of illustrating the content.  
 Images fournies uniquement à titre d'illustration du contenu.  
 الرسوم لأغراض توضيحية فقط.

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# INSTRUCTIONS FOR USE

## THESE INSTRUCTIONS ARE IMPORTANT. PLEASE KEEP THEM FOR FUTURE REFERENCE.

Dear customer,  
Thank you for choosing Pic Solution's Kappy Pic Solution aerosol therapy unit.

### INTENDED USE:

The Kappy Pic Solution piston nebuliser is meant to be used with the dedicated nebuliser cup to transform liquid drugs (in solution or suspension) into an aerosol. Thanks to accessories such as the mask, mouthpiece or nose fork the nebulised microparticles remain suspended in the air as an aerosol, which the patient then breathes in, bringing the drug directly to the patient's airways for the treatment of acute or chronic diseases according to your Physician's indications.

### INTENDED USERS

The Kappy Pic Solution can be used by adults and children of paediatric age for the aerosol treatment of chronic or acute respiratory diseases. Intended uses do not include the use of an aerosol therapy unit for the administration of lifesaving drugs.

There are no intrinsic contraindications to the use of the device: any restrictions on use depend on the type of drug used for aerosol treatment.

The device can also be used by pregnant women, for specific treatments prescribed by qualified physicians and/or healthcare personnel. It can be used by children and/or patients with reduced physical and/or cognitive abilities, so long as they receive assistance and are under adult supervision.

### GENERAL CONDITIONS FOR USE

The device is meant to be used in a household environment, it may also be used in a professional setting with the limitations listed above.

For increased safety and hygiene, never use the same accessories on more than one patient. Purchase a dedicated kit for each user, see the list of compatible accessories listed in next paragraph.

Replace the accessories after one year of use or if there are any visible signs of deterioration of the materials themselves.

The instructions for use must be read to guarantee the safe use of the device.

### ACCESSORIES

- A - Adult mask
- B - Children's mask
- C - Mouthpiece
- D - Nose fork
- E - Connection tube
- F - Nebuliser cup

Available as spare parts in kit:

**REF** 58038301000000: A, B, C, D, E, F, O

**REF** 58038302000000: A, E, F

**REF** 58038303000000: B, E, F

### DEVICE COMPONENTS

- G - On/Off switch
- H - Compressed air outlet
- I - Filter
- L - Nebuliser cup support
- M - Openings for cooling
- N - Power cord
  
- O - Filters

## GENERAL WARNINGS

- Any serious accident that occurs in relation with use of the device must be reported to the manufacturer and to the local authorities/government agencies.
- Do not use the device in conjunction with systems that provide ventilation support and are capable of releasing volatile anaesthetics or oxygen.
- The biocompatibility of the accessories and of the device's pneumatic circuit was tested in accordance with the requirements of: ISO 10993-1 - ISO 10993-5 - ISO 10993-10 ISO 10993-23 - ISO 18562-1 - ISO 18562-2 - ISO 18562-3 - ISO 18562-4.
- The masks have holes that optimise the efficiency of the aerosol treatment, do not clog them.
- Before use, carefully read the information in this manual, and keep for future reference. If there are doubts about the interpretation of any contents of this instructions manual, contact the dealer, distributor or manufacturer.
- Unit designed for aerosol therapy. Use the drugs prescribed in accordance with the instructions provided by your doctor and by the drug manufacturer. Any use other than the one for which the unit is intended is considered improper and therefore dangerous; PIKDARE cannot be held liable for any damage caused by incorrect or unreasonable use or if the device is used in electrical systems that do not conform to safety regulations in effect.
- After removing the device from its packaging, check that it is intact, without any visible damage that may have occurred during transport. If in doubt, do not use the unit, and contact the Pic Service Centre.
- Packaging (bags, box, etc.) may be dangerous and should be kept out of reach of children. Keep the device out of reach of children, pets and persons with cognitive limitations.
- Before connecting the device, check that your mains voltage corresponds to that shown on the rating label. The rating label is located on the bottom of the unit.
- In case of incompatibility between the unit plug and the electrical socket, use certified adaptors, provided applicable legislation in the individual country allows such use, or have a PIKDARE S.p.A. authorised technician or personnel from the Pic service centre replace the plug.
- To prevent overheating and damage to the compressor, switch off the appliance for at least 40 minutes after every 20 minutes of uninterrupted use.
- **WARNING! Using the device before 40 minutes have elapsed may result in motor overheating and cause the safety cut-off devices to activate.**
- To ensure correct operation, make sure that the air filter is dry. The filter cover may be removed using a tool with a flat blade. **WARNING! Keep the filter cover and the filter out of the reach of children as these are small parts and could represent a swallowing and suffocating hazard.**
- Do not use the device in the presence of nitrous oxide, oxygen or flammable air/anaesthetic mixtures.
- Turn off the device and disconnect from the mains after use and before adding more medicine. Do not fill the nebuliser cup beyond the maximum level indicated on it and check that it is assembled with all its components.
- Keep the device and power cord away from hot surfaces.
- Do not use the appliance while taking a bath or shower, or when in a damp place, or close to bath-tubs, sinks, washbasins, or in any other situation in which there are liquids that may come into contact with the appliance.
- Never touch the device with wet or damp hands.
- Do not drop or put the device into water or other liquids. Should this occur, unplug it immediately, stop using the appliance and contact PIKDARE S.p.A. or a Pic Service Centre.
- Do not block the openings for cooling air.
- Do not use the device if you feel drowsy.
- The device must not be used for patients who are unconscious or who do not breathe spontaneously.
- This product is not a toy. This product must be assembled, maintained and handled by an adult.
- If children or persons with physical or mental disabilities use the device, proper supervision by an adult is required. Children must be kept under adult supervision to ensure that they do not play with the device. **The device contains small parts that can be swallowed and the power cord poses risk of strangling.**
- **Even the tube represents a strangling hazard.**
- Do not use connectors or accessories not recommended by the manufacturer.
- Make sure the unit is placed on a level, stable surface when used in order to prevent spillage.
- During use, the surface must be free of any objects that may obstruct the correct flow of air.
- Do not pour liquid on the device, especially on the filter cover: this could lead to its undesired discharge/nebulisation.
- Before cleaning or performing any maintenance work on the device, disconnect it from the mains by pulling the

plug out of the electrical socket.

- If you decide to no longer use the device, disconnect from the mains and cut off the power cord so the device cannot work. Dispose of the device and power cord in compliance with existing laws. We also recommend eliminating any parts that are potential hazards, particularly for children.
- Refer to laws in force for information regarding the disposal of accessories subject to wear. To dispose of the device, refer to Directive 2012/19/EU.
- Do not use the device if, after a fall, it shows signs of damage on any of its parts. If in doubt, do not use the device, and contact the Pic Service Centre.
- In the case of failure or malfunction of the device, switch it off by removing the plug from the mains and contact the Pic Service Centre only. The power cord and fuse listed in the "Technical Specifications" must be replaced by Pic Service Centre staff only.
- The patient and operator are the same person. The only operations that the patient/user may carry out safely are those listed in the paragraphs "Preparation and use of the device", "Compressor cleaning and maintenance", "Cleaning and maintenance of the accessories".
- **WARNING! The power cable is the separating element from the electrical mains. When in use, to disconnect the appliance from the mains, switch it off and pull out the plug from the power outlet. Always connect the plug to an easy to reach socket to ensure that the device may be easily disconnected from the electrical mains. Do not attempt to open, repair or tamper with the appliance.**
- This activity may expose the user and the patient to an electrical hazard if not carried out by authorised personnel.
- Replace the accessories after one year of use or if there are any visible signs of deterioration of the materials they are made of.

## PREPARING AND USING THE UNIT

- Take the device and the accessories out of the package.
- Pour the medicine and saline solution in the nebuliser tank, according to the doses recommended by your physician, without exceeding the maximum level shown by the mark on the tank; close the nebuliser cup (fig. 4).  
**Note:** the notches on the nebuliser cup are for illustrative purposes only and do not provide the exact dosage of the drug contained in the cup.
- Connect the nebuliser cup to the accessory needed for the therapy: Adult or children's mask, fork, mouthpiece (applied parts, see fig. 3). In case of doubt, ask your physician about which accessory to use.
- Connect the tube to the compressed air outlet and to the nebuliser cup;
- Connect the device to the mains;
- Turn on by pressing the on/off switch and start therapy;
- During the treatment session, the nebuliser cup can be rested on its own support;
- Carry out the treatment from a relaxed, seated position keeping the nebuliser cup vertical.
- When the aerosol jet becomes intermittent, interrupt the treatment for a few seconds to allow the suspended drops of solution to fall on the nebuliser cup walls. Resume the treatment session, which will be finished when no more nebulised solution is released from the nebuliser cup;
- Once the treatment session is completed, switch the device off, unplug the device from the electrical mains, remove the accessory used, and clean the device and accessories according to the instructions in the Cleaning and Maintenance paragraph;
- Warning! Do not service and/or repair the device while it is being used or is attached to the electrical mains! Any other maintenance/service activity/operation must be done solely by a PIC Service Centre.

## COMPRESSOR CLEANING AND MAINTENANCE

Warning! Do not wet the device or put it into water or other liquids. For cleaning, use a clean, dry cloth only.

## CLEANING AND MAINTENANCE OF THE NEBULISER CUP AND ACCESSORIES

- Step 1. Remove the used accessory from the nebuliser cup.
- Step 2. Disconnect the connection tube from the nebuliser cup.
- Step 3. Disassemble all the components of the nebuliser cup (fig. 4b).
- Step 4. Discard the remaining medication.

Step 5. Wash the parts in water.

Step 6. Hand dry or air dry in a clean environment using a soft, clean, lint-free cloth.

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**WARNING! Do not boil the masks.**

## DISINFECTION PROCEDURE

Disinfect the nebuliser cup, the mouthpiece, the masks or the nose fork after the last treatment of the day as follows. Use a solution of 70% ethyl alcohol:

Step 1. Immerse the components in the disinfectant solution for 10'.

Step 2. Remove the parts and discard the solution.

Step 3. Rinse the parts with clean drinking water (60°), shake off excess water and allow to air dry in a clean environment.

Store the device and the accessories in a cool, dry place, away from light and heat.

Never use benzene, thinners or other inflammable chemical substances for cleaning. This procedure has been tested for 365 cleaning.

Clean and check that the condition of the accessories before and after use. Replace them if damaged.

## CHECKING AND REPLACING FILTERS

Periodically check the condition of the filter. The filter was inserted to protect the compressor.

Correct filter maintenance prolongs the life of the aerosol. The frequency with which filter should be replaced depends on the conditions in which the device is used. The filter is located in the upper part of the aerosol. If the air filter has changed color or has been used for more than 60 days, replace it with a new one.

### TO REPLACE THE FILTER:

- Remove the filter cover by lifting it up carefully, using a tool with a flat blade;
- Remove the filter to be replaced from its seat;
- Carefully fit the new filter;
- Reposition the filter cover by pushing it down fully in its seat, making sure that it cannot be removed without the use of a tool.

**WARNING! Never leave the removed filter cover and/or filters unattended: These are small parts posing the risk of suffocation if swallowed by children.**

| TYPE OF DEFECT      | CAUSE  | TROUBLESHOOTING  |
|---------------------|--|--|
| Scarce nebulisation | Clogged nebuliser cup                            | Remove and detach all the nebuliser cup components and clean them as illustrated in this instruction manual.                                 |
| Scarce nebulisation | Clogged nebuliser cup                            | If the issue persists after cleaning, replace the nebuliser cup.   |
| No nebulisation     | Nebuliser cup components not assembled correctly | Disassemble and reassemble the nebuliser cup; if the issue persists, replace the nebuliser cup.  |
| No nebulisation     | The compressor does not work                     | Verify that the power cord plug is properly connected and that the on/off switch is on. If the issue persists, contact a Pic Service Centre. |

|  |   |   |
|--|---|---|
| Slow nebulisation  | Medicine too thick                                      | Dilute the drug as instructed by the manufacturer.                    |
| No nebulization or low nebulization rate when the power is on. | Too much or too little medication in the nebuliser cup. | Add the correct amount of prescribed medication to the nebuliser cup. |

### TECHNICAL SPECIFICATIONS

|                                  |                                |
|----------------------------------|--------------------------------|
| Nominal voltage                  | 230V~                          |
| Frequency                        | 50 Hz                          |
| Power consumption                | 50W                            |
| Working condition                | 20'ON-40'OFF                   |
| Fuse:                            | T2AH250V                       |
| Weight of the device:            | approx. 1.2 Kg                 |
| Device service life when in use: | 3,000 cycles/treatments        |
| Service life of accessories:     | 365 cycles/treatments - 1 year |
| Noise emissions at 100 cm*       | max 58 dBA                     |
| Noise emissions at 50 cm*        | max 60.5 dBA                   |
| Max pressure:                    | 2.57 ± 0.11 bar                |
| Operational pressure:            | 0.98 ± 0.06 bar                |
| Filling (min-max):               | 2 mL-8 mL                      |

Always fill the nebuliser cup in compliance with the range defined above.

\*Data on noise levels were measured on a new device. Values may vary over time with use.

Performance:

Measurements taken with a power supply of 230V 50Hz:

|                              | <b>0.1% ALBUTEROL</b> | <b>FLUBRON SOLUTION (15 MG/2 ML OF AMBROXOL HCL)</b> | <b>PULMAXAN SUSPENSION (0.5 MG/ML OF BUDESONIDE)</b> |
|------------------------------|-----------------------|--|--|
| MMAD                         | 3,59 ± 0,18 µm        | 3,85 ± 0,06 µm                                       | 4,55 ± 0,16 µm                                       |
| GSD                          | 2,39 ± 0,11           | 2,26 ± 0,04  | 1,72 ± 0,03  |
| Respirable fraction (%<5 µm) | 61,2 ± 0,02           | 60,3 ± 0,02  | 50,3 ± 0,03  |
| % >5 µm                      | 33,3 ± 0,01           | 36,1 ± 0,01  | 38,1 ± 0,03  |
| % 2 µm to 5 µm               | 37,4 ± 0,01           | 39,8 ± 0,02  | 44,3 ± 0,01  |
| % <2 µm                      | 23,8 ± 0,03           | 20,4 ± 0,01  | 5,95 ± 0,01  |
| Output (mL)                  | 0,29 mL ± 3,52%       | 0,34 mL ± 5,15%                                      | 0,29 mL ± 4,52%                                      |
| Output rate                  | 0,10 mL ± 5,31%       | 0,11 mL ± 1,48%                                      | 0,06 mL ± 1,47%                                      |
| Residual Volume              | 1,23 mL               | 1,30 mL  | 1,22 mL  |

The data declared refers to the drugs tested and may vary depending on the drugs used, especially if these are high viscosity drugs.

The performance data listed in the table was measured using a breathing pattern that simulates adult breathing, it may vary with a child's breathing pattern.

EN ISO 27427 Respiratory therapy unit - Part 1: Nebulising systems and their components

Normal working conditions:

Temperatures: +5°C ~ +40°C

Relative humidity: 15%~85%,

Atmospheric pressure: 700 ~ 1060 hPa

Storage & transportation conditions:

Temperatures: -25°C ~ 70°C,  
Relative humidity: 15%~85%,  
Atmospheric pressure: 700 ~ 1060 hPa

\* Data on noise levels were measured on a new device. Values may vary over time with use.

This section contains specific information about product compliance with the EN 60601-1-2 standard. Kappy Pic Solution is a medical electrical device that requires special precautions regarding electromagnetic compatibility and needs to be installed and commissioned according to the electromagnetic information provided.

Mobile and portable RF communications equipment (mobile phones, transceivers, etc.) may affect the device.

Kappy Pic Solution [REF] 58038221000500 is intended for use in the electromagnetic environment specified below. The person who purchases or who uses the Kappy Pic Solution [REF] 58038221000500 is responsible for ensuring that the device is used in an environment that complies with said specifications.

#### Guidance and manufacturer's declaration - Electromagnetic emissions

| Phenomenon                                       | Professional healthcare facilities <sup>a)</sup> | DOMESTIC ENVIRONMENT <sup>a)</sup> |
|--|--|------------------------------------|
| Conducted and radiated radio frequency EMISSIONS | <sup>a)</sup>                                    | CISPR 11<br>Group 1 Class B        |
| Harmonic distortion                              | IEC 61000-3-2 <sup>b)</sup> Class A              |                                    |
| Voltage fluctuations/flicker                     | COMPLIANT WITH IEC 61000-3-3 <sup>b)</sup>       |                                    |

<sup>a)</sup> The device is intended for use in a home environment or professional healthcare structure and may exclusively be used in hospital rooms and in sections of hospitals or clinics set up for respiratory therapy. The more restrictive acceptance limits of Group 1 Class B (CISPR 11) have been considered and applied. The device is suited for use in the aforementioned environments so long as it is connected to the public mains.

<sup>b)</sup> The test may be applied to this environment provided that the ME EQUIPMENT and the ME SYSTEM used are connected to the PUBLIC MAINS and that the power supply is compliant with the prescriptions set by the reference policy governing electromagnetic compatibility (EMC).

#### Guidance and manufacturer's declaration – Electromagnetic immunity – Enclosure port

| Phenomenon  | Basic EMC standard or test method | Immunity test levels   |   |
|---|-----------------------------------|--|---|
|   |                                   | Professional healthcare facility environment   | Home healthcare environment                                 |
| ELECTROSTATIC DISCHARGE                               | IEC 61000-4-2                     | ± 8kV contact<br>± 2 kV, ±4kV ±, ±8 kV, ±15 kV air   |   |
| Radiated, radio-frequency electromagnetic fields      | IEC 61000-4-3                     | <sup>a)</sup>  | 10 V/m <sup>b)</sup><br>80 MHz - 2.7 GHz 80%<br>AM at 1 kHz |
| Fields located near RF wireless communication devices | IEC 61000-4-3                     | <b>COMPLIANT</b><br>NOTE: Further information about the distances to be maintained between portable and mobile RF communications device (transmitters) and the Kappy Pic Solution [REF] 58038221000500 can be requested from PIKDARE S.p.A. using the contact information provided in this manual. It is nevertheless advisable to keep the electromechanical aerosol therapy device at a safe distance (at least 0.5 m) from mobile phones or other communication devices that use radio frequencies, to minimise interference. |   |
| RATED power frequency magnetic fields.                | IEC 61000-4-8                     | 30 A/m <sup>c)</sup><br>50 Hz or 60 Hz   |   |

) The device is destined for household use or use in professional healthcare facilities and may only be used in the patient's room and in locations designated for respiratory therapy in a hospital or clinic. The most restrictive limits were applied for the immunity tests.

) Before modulation.

) This test level presupposes a minimum distance of at least 15 cm between the Kappy Pic Solution **REF** 58038221000500 and sources of magnetic fields at grid frequency.

#### Guidelines and Manufacturer's Declaration - Electromagnetic Immunity - Alternating Current Entry Port

| Phenomenon  | Basic EMC standard or test method | Immunity test levels  |                             |
|---|-----------------------------------|---|-----------------------------|
|   |                                   | Professional healthcare facility environment  | Home healthcare environment |
| Electrical fast transient/ burst  | IEC 61000-4-4                     | ± 2 kV<br>100 kHz repetition frequency  |                             |
| Line-to-Line surges   | IEC 61000-4-5                     | ± 0.5 kV, ±1 kV   |                             |
| Line-to-earth surges  | IEC 61000-4-5                     | Class II device or equipped with Class II AC/DC adapter. Test not applicable.<br>± 0.5 kV, ±1kV; ±2kV |                             |
| Conducted disturbances induced by radio frequency fields.   | IEC 61000-4-6                     | 3 V *) 0.15 MHz - 80 MHz<br>6 V *) for ISM bands between 0.15 MHz and 80 MHz *)<br>80% AM at 1 kHz    |                             |
| Voltage dips.   | IEC 61000-4-11                    | 0% in UT; 0.5 cycle<br>At 0°,45°,90°,135°,180°, 225°, 270° and 315°.                                  |                             |
| Voltage interruption  | IEC 61000-4-11                    | 0% in UT; 1 cycle and 70% in UT; 25/30 cycles<br>Single phase: at 0°                                  |                             |
| *) Effective value (rms) before modulation.   |                                   |   |                             |
| *) The frequency bands from 0.15 MHz to 80 MHz used by industrial, scientific and medical (ISM) radiofrequency devices are: 6.765 MHz - 6.795 MHz, 13.553 MHz - 13.567 MHz, 26.957 MHz - 27.283 MHz and 40.66 MHz - 40.70 MHz. The radio bands from 0.15 MHz to 80 MHz used for amateur radio are: 1.8 MHz - 2.0 MHz, 3.5 MHz - 4.0 MHz, 5.3 MHz - 5.4 MHz, 7 MHz - 7.3 MHz, 10.1 MHz - 10.15 MHz, 14 MHz - 14.2 MHz, 18.07 MHz - 18.17 MHz, 21.0 MHz - 21.4 MHz, 24.89 MHz - 24.99 MHz, 28.0 MHz - 29.7 MHz and 50.0 MHz - 54.0 MHz. |                                   |   |                             |

More information on distances to be kept between portable and mobile radio communications equipment and the Kappy Pic Solution **REF** 58038221000500 can be requested from PIKDARE S.p.A., through the contact information provided in this manual. It is nevertheless advisable to keep a safe distance (at least 1 m) between the aerosol therapy unit and mobile phones or RF transmitters in order to reduce possible interference.

#### LEGEND OF SYMBOLS:

|   |                                 |   |                        |
|---|---------------------------------|---|------------------------|
|  | Warning!                        |  | Humidity limitation    |
|  | Follow the Instructions for Use |  | Temperature limitation |
|  | Class II device                 |  | Recyclable packaging   |
|  | On - Off                        |  | Der Grüne Punkt        |

|   |  |   |  |
|---|--|---|--|
| <b>Hz</b>   | Frequency                              |  | Legal Manufacturer   |
| <b>SN</b>   | Serial number                          | <b>LOT</b>  | Lot Number   |
| <b>REF</b>  | Product identification code            |  | Recycling (cardboard)  |
|  | Device with Type BF applied part       |  | Triman symbol, symbol for disposal in accordance with French law   |
| <b>CE</b><br><b>1936</b>  | Complies with Regulation (EU) 2017/745 |  | Device tested by Università di Parma-Dipartimento di Farmacia to measure MMAD, aerosol output and aerosol output rate and breathable fraction. |
|  | Alternate current                      |   |  |
| <b>MD</b>   | Medical Device                         | <b>IP21</b>   | Device protected against penetration of solids and liquids. Device protected against vertically dripping water.                                |
|  | Atmospheric pressure limitation        | <b>UDI</b>  | Unique Device Identifier   |

**Device date of production:** the first pair of digits in the lot number identifies the year, the second pair identifies the month. (**LOT** 24100000; 24 = 2024; 10 = October)

#### **THIS PRODUCT COMPLIES WITH THE DIRECTIVE 2012/19/EU.**

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the unwanted appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

#### **WARRANTY**

The product is guaranteed against any lack of conformity under normal use as specified in the instructions. The warranty shall not therefore apply in the case of damage caused by improper use, wear or accidental events. For the duration of the warranty on conformity defects please refer to the specific provisions of applicable national laws in the country of purchase, where appropriate.

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