# VEDOECO PLUS

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### Intend use

Used to measure human body temperature

**Data**  
Measuring range: 35.5°C-42°C - Resolution: 0.1°C - Accuracy: +0.1°C/-0.15°C - Operating temperature: -10°C to 35°C - Storage temperature: -10°C to 42°C

### Quality guarantee

Mercury-free medical liquid-in-glass thermometers are clinical thermometers with maximum device. They are intended for measuring body temperature. The manufacture and testing of the thermometers is conducted in compliance with the standards: -EU council directive 93/42/EEC on Medical devices. - EN 12470-1+A1:2009 Clinical thermometers

-Part 1, Metallic liquid-in-glass thermometers with maximum device. A quality management system certified under EN ISO 13485 guarantees compliance with the required quality criteria.

### Precautionary measures

1. Do not drop the thermometer or subject it to sudden impacts.
2. Do not bend or bite the tip of the thermometer.
3. Keep the thermometer out of children's reach.

4. Thermometers with a damaged outer casing must not be used (risk of injury).
5. Avoid high temperature storage, not more than 42°C.

- Before measuring temperature, you should avoid drinking hot or cold drinks, drinking alcohol, or strenuous exercise that may affect your body temperature.

- Before use, check the product integrity (absence of liquid column disconnection and crack), verify the liquid column stability and paper position (the small line on the glass should be aligned with the 38°C indicator). If damaged, do not use the product.

- Do not use the thermometer if the liquid column drops significantly within 1 minute after its removal from the measuring site. Please contact the distributor or the Italian toll-free number 800 900 080 for after-sales.

- In case of any accident during use, call the emergency number immediately and seek the help of medical professionals.

- Children should be supervised by an adult when using it to prevent biting and swallowing.
- When using it, use it best to put the product temperature bubble downward in the original case.

- Do not use over a wound.

### Instructions for use

1. Before taking the patient's temperature, check whether the column of liquid is below the first numbered line (between 35°C and 36°C).

2. Clean the glass thermometer before and after each use. In doing so, please pay attention to the "Note on Cleaning" below.

3. Position the thermometer on the body according to the method of measurement to be used (oral, rectal or axillary).

4. After about 5 min the thermometer will have reached the body temperature. The temperature can be read by rotating the thermometer gently. The measuring fluid in the glass capillary of the thermometer shows the measured body temperature on the temperature scale.

5. After measuring, the liquid must return below 36°C. To do so, put the thermometer back into its casing and make sure to push it to the bottom by passing the probe through the relevant hole (Fig. 1). Then take the casing and shake it with swift, rhythmic movements to make the liquid fall below the first numbered notch.

- **NOTE:** Always disinfect the thermometer before use. In hospitals disinfection must be carried out by experienced personnel.

### Professional methods of temperature measurement

#### Oral

Oral measurement is preferred by physicians. It provides accurate results quickly. The thermometer bulb (probe) should be placed in one of the two thermal pockets to the left or right of the root of the tongue. The normal temperature using this measurement method is 37 °C (measuring time 5-7 minutes).

#### Rectal

Rectal measurement is the method of choice for use with young children. The thermometer bulb (probe) is inserted into the anus.

**NOTE:** Do not insert the tip of the probe more than 2 cm into the rectal; stop if you meet with any resistance. A little skin cream or baby oil may be applied to the tip if necessary. Do not use Vaseline (petroleum jelly). The normal temperature using this method is approximately 37.6°C (measuring time 5-7 minutes).

#### Ampit

Rub the amplit with a dry cloth. The thermometer vesel should now be placed in the amplit and the arm pressed against the side of the body. The normal temperature here is 36,7°C (measuring time 7-10 minutes).

**NOTE**  
1. In children, an adult should be present throughout the whole process of measuring and recording the temperature.  
2. When measuring, please take the necessary care to avoid any danger of injury.  
3. All thermometers have been tested by the manufacturer to ensure that they comply with the prescribed tolerance threshold of (+0.1/-0.15°C).

### NOTES ON CLEANING

The thermometer has a hermetically sealed, hygienic glass casing and can be disinfected without restriction. The temperature of the cleaning and disinfectant solution must not exceed 42°C. The cleaning of the glass thermometer before and after use can be done with cold soapy water. Disinfect the thermometer using a piece of cotton wool or cloth that has been soaked in alcohol or dip the thermometer in a disinfectant solution. In hospitals, disinfection should always be done by expert personnel. If the thermometer breaks, any spilled liquid can be removed using a brush, tissue or cloth dampened with soapy water. The broken pieces may be disposed of as household waste.

### Limited warranty

Subject to normal use, this thermometer is guaranteed against any manufacturing defects for a period of two years from the date of production. Damage to the thermometer caused by improper use is not covered by the warranty.

### RE-CALIBRATION IS NOT NECESSARY!

### Notice to the user or patient

Please be noticed that, any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### SYMBOL INDEX

	Follow the instructions		Temperature limitations
	Caution		Store in a dry place
	Batch code (mm/yyyy; month/year) on the reverse of the outer casing (applied by means of laser without LOT symbol)		Fragile: handle with care
	Manufacturer		Triman icon
	Production date		Der Grüne Punkt
	Symbol of Compliance with Directive 93/42/CEE		Recyclable packaging
	Distributor		Keep away from sunlight
	Mandatario europeo		Importatore

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### CONSIGNES DE NETTOYAGE

The thermometer has a hermetically sealed, hygienic glass casing and can be disinfected without restriction. The temperature of the solution detergent and disinfectant should not exceed 42°C. The cleaning of the glass thermometer before and after use can be done with cold soapy water. Disinfect the thermometer with a piece of cotton wool or cloth that has been soaked in alcohol or dip the thermometer in a disinfectant solution. In hospitals, disinfection should always be done by expert personnel. If the thermometer breaks, any spilled liquid can be removed using a brush, tissue or cloth dampened with soapy water. The broken pieces may be disposed of as household waste.

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