# Mediblink®

medical devices

## **EN** MEDIBLINK Compressor Nebulizer Airy M430

### INSTRUCTIONS FOR USE

We encourage you to thoroughly read this guidebook to learn about the features of this product.

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# TABLE OF CONTENTS

<b>1. Introduction</b>	6
1.1 Intended use/Indication for use	6
<b>2. Product identification</b>	6
<b>3. Important safeguards</b>	7
3.1 Product warnings	7
3.2 Operating warnings	8
3.3 Storage warnings	9
3.4 Cleaning warnings	9
<b>4. Operating your compressor nebulizer</b>	9
<b>5. Cleaning</b>	10
5.1 Rinsing the accessories (after each treatment)	10
5.2 Cleaning the compressor	10
5.3 Changing the Air Filter	10
<b>6. Disinfection</b>	11
<b>7. Particulate distribution diagram</b>	11
7.1 Particulate distribution diagram	11
7.2 Particulate distribution data	12
<b>8. Specifications</b>	12
<b>9. EMC information</b>	13
<b>10. List of symbols</b>	17
<b>11. Warranty</b>	18

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# 1. Introduction

---

Thank you for purchasing the Mediblink Compressor nebulizer Airy M430. It is a compact medical device designed to efficiently deliver physician prescribed medication to the bronchial lung passages. With proper care and use, it will provide you with many years of reliable treatment.

## 1.1 Intended use/Indication for use

It is intended for use in the treatment of asthma, COPD and other respiratory ailments in which an aerosolized medication or saline solution is required during therapy. The device is for pediatric and adult population.


Mediblink Airy M430 nebulizer should be used under the supervision of a licensed physician and/or a respiratory therapist. We encourage you to thoroughly read this guidebook to learn about the features of this product. Any use of this product other than its intended use should always be avoided.

## 2. Product identification

---

Mediblink Airy M430 nebulizer comes with the following components:

- a) Main unit and parts (Figure 1)
  - Compressor
  - Air filter
  - Power cord
  - User manual
- b) Accessories (Figure 2)
  - Nebulizer cup
  - Air tube
  - Mouthpiece
  - Big/Adult mask and/or Small/Pediatric mask

**Figure 1:** Main unit and parts


Figure 2: Accessories

## 3. Important safeguards

**Note:** Please read this user manual thoroughly before using the product. The following basic precautions are needed when using an electrical product:

### **Warning**

Failure to read and observe all precautions could result in personal injury or equipment damage.

### 3.1 **Product warnings:**

- a) To avoid electrical shock:
  - Do not connect the power cord terminal by using wire or other metallic conductors.
  - Do not use the compressor (main unit) or power cord while they are wet.
  - Keep the device away from water.
  - Do not use or store the device in humid locations, such as a bathroom. Use the device within the operating temperature and humidity.
  - Use only the power cord designed by manufacturer for this device. Use of any other power cord may damage the device.
  - Do not immerse the power cord or the main unit in liquid.
  - Do not use while bathing.
  - Do not reach for a device that has fallen into water - immediately unplug the device.
  - Do not overload power outlets. Plug the power cord into the appropriate voltage outlet.
- b) Never operate the device if it has any damaged parts (including power cord), if it has been dropped or submersed in water. Promptly send it to a service center for examination and repair.
- c) The device should not be used where flammable gas, oxygen or aerosol spray products are being used.
- d) Disconnect the device from the electrical outlet before cleaning, filling and after each use.

## 3.2 Operating warnings:

- a) Connect this product to an appropriate voltage outlet for your model.
- b) Never operate it if this device has a damaged cord or plug, if it has been dropped or in any way submersed into water, or if it does not work properly. Return it to a service center for repair.
- c) If any abnormality occurs, discontinue use immediately until the device has been examined and repaired.
- d) If you feel anything unusual during use, stop using the device immediately and consult your physician.
- e) Do not cover the device with a blanket, towel, or any other type of cover during use. This could result in the compressor overheating or malfunctioning.
- f) Always unplug the device immediately after use.
- g) Never block the air openings of the main unit or place it where the air openings may be obstructed.
- h) The compressor and power cord are not waterproof. If liquid spills on these parts, immediately unplug the power cord and wipe off the liquid with gauze or other soft absorbent material.
- i) Do not use a cellular phone near the device. It may result in an operational failure.
- j) Do not use extension cords. Plug the power cord directly into the electrical outlet.
- k) For type, dose, and regime of medication follow the instructions of your licensed physician and/or respiratory therapist.
- l) Do not add more than 8 ml of medication to the medication cup.
- m) Do not use the device if the air tube is bent.
- n) Clean and disinfect the nebulizer, air tube, mouthpiece or masks before using them for the first time after purchase.
- o) Provide close supervision when this device is used by, on, or near infants, children or compromised individuals.
- p) Make sure the air filter is clean. If the air filter has changed color or has been used on average for more than 30 days, replace it with a new one.
- q) Make sure the air tube is correctly connected to the compressor and the nebulizer. Air may leak from the air tube during use if not securely connected.
- r) Use only authorized parts and accessories from the manufacturer. Parts and accessories not approved for use with the device do not perform the expected specifications or they may damage the device.
- s) To avoid the medication residue on the face after nebulization, be sure to wipe the face after removing the optional masks.
- t) Do not use the device where the device may be exposed to flammable gas or vapors.
- u) Using a solution, suspension, or emulsion different from that recommended by the manufacturer, in particular, a suspension and/or high-viscosity solution, can alter the particle size distribution curve, the mass median aerodynamic diameter (MMAD), aerosol output, and/or aerosol output rate, which can then be different from those disclosed by the manufacturer.
- v) The nebulizer and other accessories are designed for single-patient reuse, do not reuse them among different people to avoid cross-contamination.
- w) The nebulizing system is not suitable for use in an aesthetic breathing system or a ventilator breathing system.

### 3.3 ⚠ Storage warnings:

- Do not store the device under direct sunlight, high temperature or humidity. Store the device within the storage temperature and humidity.
- Keep the device out of reach of small children.
- Always keep the device unplugged while not in use.
- Store the device and the components in a clean and safe place.
- Do not wrap the power cord around the compressor (main unit) and power cord.


### 3.4 ⚠ Cleaning warnings:

- Do not immerse the device in water. It may damage the device.
- Disconnect the device from the electrical outlet before cleaning.
- Clean all necessary parts after each use as instructed in this user manual.


## 4. Operating your compressor nebulizer

**Note:** Prior to initial operation, the nebulizer should be thoroughly cleaned referring to the “Cleaning Procedures” in the manual.

- Place your nebulizer on a flat and stable surface. Be sure that you can easily reach the controls when you are seated.
- Open the clear cover and remove the inside accessories.
- Gently twist the top part of the nebulizer cup counter-clockwise to disassemble the nebulizer cup.
- Fill the bottom section of the nebulizer cup with the medication prescribed by your physician. Be sure that the cone is put inside the bottom section (Figure 3).
- Gently twist the top part clockwise to re-assemble the nebulizer cup. Be sure the two sections fit well.
- Attach one end of air tube to the base of the nebulizer cup (Figure 4).
- Attach the other end of air tube to the air tube connector located on the front of the compressor.

**Figure 3:** Nebulizer cup

**Figure 4:**  
Attaching the accessories

## **Warning**

- Make sure the nebulizer cup is correctly assembled, the air filter is properly installed, and the air tube is correctly connected to the compressor and the nebulizer cup. Air may leak from the air tube during use if not securely connected.
  - Do not use the device if the air tube is bent.
- h) Attach mouthpiece or mask per your choice to the top section of nebulizer cup.
  - i) Plug the power cord into an appropriate electrical outlet. Make sure at this stage, the power switch is at “OFF” status.
  - j) Press power switch to begin your prescribed treatment.

**Important:** The compressor motor has a thermal protector which will shut off the device before the device is overheated. When the thermal protector shuts the device off, please:

- a) Switch off the device.
  - b) Unplug the device from the electrical outlet.
  - c) Wait 30 minutes for the motor to cool down before another treatment. Make sure the air openings are not obstructed.
- k) When treatment is finished, shut off the device and unplug it from the electrical outlet.

## 5. Cleaning

---

### 5.1 Rinsing the accessories (after each treatment)

- a) Disconnect the air tube, nebulizer cup, mouthpiece or mask.
- b) Gently twist the nebulizer cup to open it.
- c) Discard remaining medication in the nebulizer cup.
- d) Rinse the nebulizer cup, mouthpiece or mask with tap water.
- e) Dry them with clean soft cloth, then let it air dry.
- f) Re-assemble the nebulizer cup when completely dry and put these parts in a dry, sealed container.

### 5.2 Cleaning the compressor

- a) Clean the casing of the main unit by using a soft cloth moistened with water or a mild detergent.
- b) Do not use abrasive cleaners.
- c) Dry the casing immediately by using a soft, dry and clean cloth.

### 5.3 Changing the air filter

- a) Do not use cotton or any other materials. Do not wash or clean the filter. Use only filters supplied by Mediblink. And do not operate without a filter.
- b) Change the filter every 30 days or when the filter turns gray.
- c) Before inserting the new air filter, make sure the air filter is clean and free of dust. Do not operate the device without the air filter. Use only the Mediblink air filter designed for this device.
- d) Changing procedure:

- i) Pull the air filter cover to remove from the back side of the compressor.
- ii) Remove the dirty air filter.
- iii) Replace the used filter with a new one.
- iv) Put the air filter cover back on the compressor.

## 6. Disinfection

---

Please follow the following steps to disinfect your accessories (including mouthpiece, mask, air tube, and nebulizer cup) unless otherwise specified by your physician. It is suggested that the accessories are disinfected after the last treatment of the day.


- a) After cleaning the accessories, submerge them into 75% medicinal alcohol solution for 10 minutes. Make sure the accessories are immersed into the liquid level completely.
- b) Then let the accessories air dry on a dry & clean cloth.

## 7. Particulate distribution diagram

---

### 7.1 Particulate distribution diagram

The following table is particulate distribution diagram, the test particles are respective NaCl.

**Figure 5:** Particulate distribution diagram of NaCl

## 7.2 Particulate distribution data

Solution: albuterol 0.1% (M/V) concentration in 0.9% sodium chloride solution.

Item		Adult	Pediatric
Percentage of particulate distribution (PORR)	>5 $\mu\text{m}$	60%	60%
	>2 $\mu\text{m}$ and <5 $\mu\text{m}$	75%	75%
	<2 $\mu\text{m}$	80%	80%
MMAD		approx. 3.65 $\mu\text{m}$	approx. 3.65 $\mu\text{m}$
Nebulization rate		0.20 ml/min	0.20 ml/min
GSD		1.798	1.798
Percentage of fill volume emitted (8ml)		85%	85%
Residual volume		1.5 ml	1.5 ml
<p><b>Note:</b> Nebulizer performance is based upon testing that utilizes adult ventilatory patterns and are likely to be different from those stated for pediatric populations.</p>			

## 8. Specifications

Power supply	AC 230V 50 Hz
Rated power	$\leq 65$ W
Air Flow	$\geq 8$ L/min
Medication Capacity	2–8 ml
Particle Size (MMAD)	approx. 3.65 $\mu\text{m}$
Percentage of particulate distribution (PORR)	approx. 80%
Nebulization rate	approx. 0.20 ml/min
Sound Level (noise)	$\leq 55$ dB
Working pressure	50–150 KPa
Water protection	IPX0
Working temperature	10 °C to 40 °C
Working humidity	30% to 75% RH
Working atmospheric pressure	700–1060 hPa
Storage Temperature	-25 °C to 70 °C
Storage conditions	10% to 95% RH
Transportation conditions	-25 °C to 55 °C, 10% to 95% RH
Standard requirements	Comply with EN 60601-1-2, EN60601-1 and EN60601-1-11

## 9. EMC information

The device needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and the device can be affected by portable and mobile RF communications equipment.

### **EMC Warnings:**

- 1) The portable nebulizer is intended to be used in the hospital and family setting except for near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.
- 2) Do not use a mobile phone or other devices that emit electromagnetic fields near the device. This may result in incorrect operation of the device.
- 3) **Warning:** This device has been thoroughly tested and inspected to assure proper performance and operation.
- 4) **Warning:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this device and the other equipment should be observed to verify that they are operating normally.
- 5) Use of accessories other than those specified or provided by the manufacturer of this nebulizer could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and in improper operation.
- 6) Do not expose the device to RFID systems.

### Essential performance of the Mediblink Compressor Nebulizer Airy M430:

Reference number	M430
Model Number	JLN-2300AS
Power on/off	Normally turn on and turn off
Air Flow	$\geq 6.5$ L/min

#### Guidance and manufacturer's declaration – electromagnetic emissions

The M430 (JLN-2300AS) is intended for use in the electromagnetic environment specified below. The customer or the user of the M430 (JLN-2300AS) should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The M430 (JLN-2300AS) uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The M430 (JLN-2300AS) is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Complies	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

### Guidance and manufacturer's declaration – electromagnetic immunity

The M430 (JLN-2300AS) is intended for use in the electromagnetic environment specified below. The customer or the user of the M430 (JLN-2300AS) should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV Contact ±15 kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. If ESD interfere with the operation of equipment, counter measurements such as wrist strap, grounding shall be considered.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for Power supply lines ±1 kV for Input/output lines	±2 kV for Power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% $U_T$ ; 0.5 cycle at 0, 45, 90, 135, 180, 225, 270 and 315 degrees, 0% $U_T$ ; 1 cycle 0 degree 70% $U_T$ (30% dip in $U_T$ ) for 25 cycle single phase at 0 degree 0% $U_T$ ; 250 cycle	0% $U_T$ ; 0.5 cycle at 0, 45, 90, 135, 180, 225, 270 and 315 degrees, 0% $U_T$ ; 1 cycle 0 degree 70% $U_T$ (30% dip in $U_T$ ) for 25 cycle single phase at 0 degree 0% $U_T$ ; 250 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of M430 (JLN-2300AS) requires continued operation during power mains interruptions, it is recommended that the M430 (JLN-2300AS) be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz/60 Hz	30 A/m 50 Hz/60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

## Guidance and manufacturer's declaration – electromagnetic immunity – For EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

### Guidance and manufacturer's declaration – electromagnetic immunity

The M430 (JLN-2300AS) is intended for use in the electromagnetic environment specified below. The customer or the user of the M430 (JLN-2300AS) should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM bands between 150 kHz to 80 MHz 80% AM at 1 kHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the M430 (JLN-2300AS), including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz	10 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2,7 GHz  where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

## Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM – For EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

### Recommended separation distances between Portable and mobile RF communications equipment and the JLN-2300AS

The M430 (JLN-2300AS) is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the M430 (JLN-2300AS) can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the M430 (JLN-2300AS) as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1.2 \sqrt{P}$	$d = 1.2 \sqrt{P}$	$d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## 10. List of symbols

Symbol	Explanation
	Electrical devices are recyclable material and should not be disposed of with household waste after use! Help us protect the environment and save resources and take this device to the appropriate collection points. Please contact the organization which is responsible for waste disposal in your area if any questions.
	Applied part of type BF
	Class II appliance
	Before use, read the instructions. Electronic instructions for use: <a href="http://www.mediblink.com/ff/m430.pdf">http://www.mediblink.com/ff/m430.pdf</a>
	Manufacturer
	Authorized Representative in EU
	Warning!
	Batch number*
	Product reference number
	Model number
	Medical device
	Unique Device Identifier
	Number of products included in packaging
	CE Marking of Conformity

### Degree of safety in the presence of flammable anesthetics or oxygen:

No AP/APG (Not suitable for use in the presence of flammable anesthetic or oxygen)

\* The production date can be read from LOT number [YYYYMM]; first four digits represent the Year and last two digits Month of production.  
Example: LOT 202503 = March 2025

# 11. Warranty

---

**Product: Mediblink Compressor nebulizer Airy M430**

**Manufactured for (Importer for EU and Distributor):** Mediblink d.o.o., Gubčeva cesta 19, 8210 Trebnje, Slovenia; info@mediblink.com; www.mediblink.com

**Seller's name, address, signature and stamp\*:**

**Date of extradition/sales\*:**

\* If the invoice is accompanied by this warranty, and if all above information can be seen from the invoice, it is not necessary to fill in this field.

## Warranty terms

Dear customers!

The warranty period is **5 years** and starts from the day of product purchase. The warranty period on spare parts: masks, mouthpiece, tube and filters is **1 year**. In case of product claim, you have to show the invoice. We kindly ask you to save the invoice!

Unfortunately, wrong handling with the device is a reason for 95% of customer complains. You can easily avoid any problem, by getting useful information provided by our special service department. To reach our service department, you can call or send e-mail to Mediblink local Distributor.

Before sending the product back to retailer, we kindly ask you to call our service department, to get help about how to use the device to save you with unneeded trips.

The Manufacturer guarantees free elimination of all imperfections due to defects in material or manufacturing procedure by repairing or replacing the product. In case that the product cannot be repaired or replaced, the customer will get the money refund. The guarantee is not valid in case of the force majeure, accidents or unexpected events (such as lightning, water, fire etc.), incorrect use or incorrect transport, non-compliance with safety and maintaining regulations or in case of unprofessional product intervention.

Traces of every day product usage (scratches, abrasions) and not subject to claim. The warranty does not eliminate the customer rights, which originate from seller responsibility for product flaws. By accepting the claimed product by the service department, the service department does not take responsibility for loss of saved data or settings on the product. All product repairs, which are performed out of product warranty period, have to be paid by customer by prior notice.

The Manufacturer guarantees the product quality and flawless product operation in the warranty period, which starts with the day of product purchase. If the product cannot be

repaired in 45 days, the product will be replaced with a new one. In case that the product cannot be replaced, the money will be refunded to the customer.

In case of product claim, call or send e-mail to Mediblink local distributor.

**Replacement parts can be ordered at [www.mediblink.com](http://www.mediblink.com).**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the Competent Authority of the Member State in which the user and/ or patient is established.

**Manufacturer:**  
Shenzhen Homed Medical  
Device Co., Ltd.  
Room 1601-1614, A Building,  
Jinxu III, No. 85 Hudi Pai,  
Songxuan Community, Guanhu  
Street, Longhua District,  
Shenzhen 518109, P. R. China

**Manufactured for**  
(Importer for EU & Distributor):  
Mediblink d.o.o.,  
Gubčeva cesta 19,  
8210 Trebnje, Slovenia;  
[info@mediblink.com](mailto:info@mediblink.com);  
[www.mediblink.com](http://www.mediblink.com)

Shanghai International Holding  
Corp. GmbH (Europe),  
Eiffestrasse 80, 20537 Hamburg,  
Germany

Instructions for use,  
version No: Ver. 1.0  
Issue date: 20. 10. 2024  
Date of last change:  
25. 10. 2024

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