## BATTERY REPLACEMENT

The battery is empty and needs replacing when the battery symbol appears on the right of the LCD. Remove the battery cover and replace it with a battery (preferably non-mercury) of the same type.

**Please note:** the «+» sign up and «-» sign down. We advise you to remove the batteries if the device is not going to be used for a longer period of time.

## PRODUCT DISPOSAL

Please ensure environmental protection. Batteries do not belong in the domestic waste. Please hand them in at collection point or the municipal recycle material centre as special waste. The alkaline battery or fuel cell may lead to excessive temperatures, fire or explosion.

This symbol on products and/or accompanying documents means that consumed electronic products must not be mixed with conventional domestic waste. Take these products to the corresponding collection points for correct treatment and recycling, where they will be accepted free of charge. For more information on the closest collection point, please enquire with your local authorities.

## TECHNICAL DATA

**Type:** Maximum thermometer.  
**Measurement range:** (32.0~42.9)°C/(89.6~109.2)°F  
**Measurement accuracy:** ±0.1°C/0.2°F (35.5°C~42.0°C/95.9°F~107.6°F).  
±0.2°C/0.4°F (for the remaining measurement range).

**Storage/transportation temperature:** (-25~55)°C, ≤95% RH.

**Ambient temperature during use:** (5~35)°C, ≤80% RH.

**Min Scale:** 0.1°C/0.1°F.

**Battery type:** Alkaline manganese battery, type LR41, 1.5 V, service life minimum 100 hours under continuous operation.

**Weight:** approx. 13 g.

## EXPLANATION OF SYMBOLS

	Battery is empty
	Product disposal instructions for electronic devices
	Manufactured by
	Manufacturing date
	Store dry
<b>IP22</b>	Classification according to the degree of protection against ingress of water as detailed in IEC 529
	Protect from direct sunlight
	The battery in this product complies with the requirements stated in European Directives 2006/66/EEC
	Type BF equipment
	LOT number
	Caution
	Follow instructions for use
<b>Lo °C or Lo °F:</b>	Temperature under 32 °C or 89.6 °F
<b>Hi °C or Hi °F:</b>	Temperature over 42.9 °C or 109.2 °F
<b>CE</b>	CE conformity marking

## LEGAL REQUIREMENTS AND GUIDELINES

This product complies with the European Directive for Medical Device 93/42/EEC and carries the CE mark, the device also complies with the specifications of below standard for:

- EN 12470-3: Clinical thermometers. Performance of compact electrical thermometers (non-predictive and predictive) with maximum device.
- EN 60601-1.
- EN 60601-1-11.
- EN 60601-1-2.

The CE Marking confirms that this is a medical device with a measuring function in the sense of the medical device act which has undergone a conformity assessment procedure. A Notified body confirms that this product fulfills all the appropriate statutory regulations.

## CALIBRATION CHECK

This thermometer is initially calibrated at the time of manufacture. If this thermometer is used according to the operation instruction, periodic re-adjustment is not required. The calibration check has to be carried out immediately, if there are indications that the product does not keep the defined error limits or the calibration properties could have been affected by an intervention or by any other means. Please also observe any national statutory regulations. The calibration check can be carried out by the competent authorities or by authorised service providers. A test instruction for calibration check can be provided to the relevant authorities and authorized services providers on request.

## WARRANTY

This product has a warranty of 2 years from the date of purchase. Please retain your receipt for the potential warranty claim. Damage resulting from incorrect use or abuse is not covered by the warranty. Battery and packaging are excluded from the warranty. Claims beyond this, including claims for damages, are excluded. If you find that the thermometer is defective and not in good function, please firstly check the battery before sending in for repair. If necessary, contact the point of sale where the thermometer was purchased.

Hangzhou Hua'an Medical & Health Instruments Co., LTD.  
Building 2, 1# Fuzhu Nan RD, Wuchang Town,  
Yuhang District 310023, Hangzhou, Zhejiang, China

Shanghai International Holding Corp. GmbH;  
Eiffelstrasse 80, 20537 Hamburg, Germany

EU: PXG Pharma GmbH,  
Pflingstweidstrasse 10-12, 68199 Mannheim, Germany.  
www.pxgpharma.com

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